

Rapid and Robust Protein Chromatography

Chromolith[®] WP 300 RP-18 2 mm I.D. HPLC Columns

The most hydrophobic of the Chromolith[®] WP 300 line, the RP-18 column is useful for the resolution of peptides and smaller proteins. One critical quality attribute (CQA) required by regulatory bodies is the peptide map of a biotherapeutic. Peptide maps generated by RP-HPLC provide valuable information about protein structure, stability, and purity. To be effective, the RP-HPLC column must be able to resolve a high percentage of the peptides in the sample. The more peptides, the better the information. The Chromolith[®] WP 300 RP-18 column gives unsurpassed RP-HPLC resolution of peptide maps from enzymatic digests. The improvements in silica and bonded-phase chemistry incorporated into the Chromolith[®] WP 300 line improve resolution by increasing efficiency and reducing peak tailing.

Key Benefits:



- Monolithic skeleton permits high flow rates to maximize throughput
- 300 Å mesopores permit large molecules to enter without fear of size-exclusion effects
- Matrix-loaded samples can be injected onto the column with little to no sample prep.

Protein Analysis on Chromolith[®] WP 300 RP-18

For large molecule separations, high efficiency separations are necessary in order to achieve resolution and good peak shape. Moving to sub-2 µm FPP-packed columns or 2.0 µm SPP-packed columns can deliver that desire; however, this comes at the cost of elevated backpressure. Chromolith[®] WP 300 RP-18, 2 mm I.D. columns provide UHPLC efficiencies, but at nearly 1/10th the backpressure.

Chromatographic conditions:

| | |
|----------------------------|--|
| Columns: | Chromolith [®] WP 300 RP-18 100-2 mm (1.52370.0001) SPP, C18, 160 Å, 2.0 µm, 100-2.1 mm FPP, C18, 300 Å, 1.7 µm, 100-2.1 mm |
| Mobile phase: | A: water (0.1% TFA) B: acetonitrile (0.08% TFA) |
| Gradient: | 4% B to 60% B in 10 minutes |
| Flow rate: | 0.38 mL/min |
| Detection: | UV, 220 nm |
| Column temperature: | 30 °C |
| Injection volume: | 0.5 µL |
| Sample: | HPLC Protein Mix 1 mg/mL, water 1) Ribonuclease 2) Cytochrome C 3) Holo-Transferrin 4) Apomyoglobin |




Excellent Lot-to-Lot Reproducibility

Chromolith® WP 300 RP-18 columns exhibit excellent batch-to-batch reproducibility, as demonstrated below with the same peptide map generated for Cytochrome C using three different Chromolith® WP 300 RP-18 columns across three different batches.

Chromatographic conditions:

| | | | |
|--------------------------|--|-----------|-----------|
| Column | Chromolith® WP 300 RP-18 100-2 mm (1.52370.0001) | | |
| Mobile Phase | A: acetonitrile 0.08% (v/v) TFA B: water 0.1% (v/v) TFA | | |
| Gradient: | Time | %A | %B |
| | 0 | 5 | 95 |
| | 25.0 | 30 | 70 |
| | 30.0 | 30 | 70 |
| Flow rate: | 0.190 mL/min | | |
| Pressure: | 18 bar | | |
| Detection: | Vanquish DAD 20 Hz, UV, 214 nm | | |
| Detector cell: | LightPipe 10 mm | | |
| Temperature: | 30 °C | | |
| Injection volume: | 0.2 µL | | |
| Sample: | Rapid Trypsin Digestion with SOLu-Trypsin Rapid Digestion Kit 2.5 mg Cytochrome C was added in a PCR vial and dissolved in 320 µL Rapid Trypsin Digestion Buffer. In the solution was 80 µL SOLu-Trypsin added and incubated at 60 °C for 1 hour in a Thermomixer. The digestion was quenched by adding 12 µL of hydrochloric acid 32 %. | | |




Antibody Fragment Analysis

Fragment analysis of a monoclonal antibody (mAb), also called middle-up analysis, is a useful technique in characterizing mAb domains without the inherent complexity of a peptide map. High efficiency is needed here to resolve subtle, structural variants of the mAb domains. The Chromolith® WP 300 RP-18 column is able to achieve the same separation efficiency and sensitivity as sub-2 µm FPP and 2.0 µm SPP-packed columns but at only 20% of the backpressure of those columns.

Chromatographic conditions:

| | | |
|--------------------------|---|-----------|
| Column | Chromolith® WP 300 RP-18, 2 mm I.D. (1.52370.0001) SPP, C18, 160 Å, 2.0 µm, 100-2.1 mm FPP, C18, 300 Å, 1.7 µm, 100-2.1 mm | |
| Mobile Phase | A: Water (0.1% (v/v) TFA) B: Acetonitrile (0.08% (v/v) TFA) | |
| Gradient: | Time (Min) | %B |
| | 0 | 20 |
| | 1 | 20 |
| | 9 | 45 |
| Flow rate: | 380 µL/min | |
| Detection: | UV, 220 nm | |
| Temperature: | 80 °C | |
| Injection volume: | 1.0 µL | |
| Sample: | SigmaMAb, 2 mg/mL (SiLu™ Lite Universal Antibody) | |
| DTT digest: | 60 µL of 40 mM Dithiothreitol (DTT) solution was added in a PCR vial, 40 µL mAb was added and incubated at 37 °C for 30 minutes creating light chain (LC) and heavy chain (HC) parts of the antibody. | |



Ordering Information

| Part Number | Description | Length (mm) | I.D. (mm) |
|--------------|--|-------------|-----------|
| 1.52370.0001 | Chromolith® WP 300 RP-18 Column | 100 | 2 |
| 1.52371.0001 | Chromolith® WP 300 RP-18 Column | 50 | 2 |
| 1.52372.0001 | Chromolith® WP 300 RP-18 Guard Columns (3 units) | 5 | 2 |

To place an order or receive technical assistance

Order/Customer Service: [SigmaAldrich.com/order](https://www.sigmaaldrich.com/order)
 Technical Service: [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice)
 Safety-related Information: [SigmaAldrich.com/safetycenter](https://www.sigmaaldrich.com/safetycenter)

Merck KGaA
 Frankfurter Strasse 250
 64293 Darmstadt, Germany

[SigmaAldrich.com](https://www.sigmaaldrich.com)

