

## Validating Cellesce's bioreactor technology for the expansion of patient-derived breast cancer organoids

Cellesce develops bioprocessing systems for the expansion of organoids in significant volumes in a way that **minimises** handling time and **maximises** reproducibility, with the long-term goal of positioning organoid technology as a cost-effective model for early-stage drug discovery.

Following the successful expansion and commercialisation of 10 human colorectal cancer organoid lines using proprietary bioreactor technology, Cellesce is now broadening its product range to include other tissue types. Through an **InnovateUK** funded project, and in collaboration with **Cardiff University**, Cellesce has been working to optimise its patented expansion process to allow the large scale expansion of **patient-derived breast cancer organoids**.

### Experimental set up

An ER-/PR+/Her2- patient-derived breast cancer organoid line (BrC-001) derived externally was used for this study. Two further PDX organoid lines (BrC-002 and BrC-003) were also derived in-house from supplied PDX tumour material, and used as controls where indicated.


All organoids were initially cultured manually according to previously described protocols (Sachs et al., 2018).

Here, in a proof of principle study, Cellesce has demonstrated the stability of breast cancer organoid phenotype, genotype and drug response following bioreactor expansion. Moreover, bioreactor-expanded organoids demonstrate reduced well-to-well variability and greater batch-to-batch reproducibility than their manually grown counterparts in drug response assays, thus reinforcing their utility in high throughput drug screening.

Cellesce will use the results of this study to direct the expansion of further breast cancer organoid lines in sufficient quantities to produce frozen batches that can be provided for sale for use 'off the shelf'.

To compare outputs from the two culture methods, organoids were divided between **(a)** manually expanded culture conditions and **(b)** Cellesce's patented bioreactor expansion process (**Figure 1**).

Organoids were cultured for the same duration in each process before harvesting from their matrices by the same protocols.








Figure 4. Gene expression analysis of organoids using RNAseq.  
Heatmap of FPKM reads of PAM50 genes in all samples. The colour scale is set by row (gene) minimum and maximum FPKM reads, where red is a higher number and blue is lower. Bioreactor samples are labelled in green for distinction.

Drug treatment assays performed on BrC-001 organoids expanded manually and using bioreactor technology, according to the workflow shown in **Figure 5**, broadly indicated that bioreactor expansion does not alter drug response; in fact in the majority of cases, analysis indicated that individual drug EC50s were not statistically different

between the two organoid sets (**Figure 6A**, indicated by \*).

Moreover, dose response data indicated that bioprocessed organoids generated tighter error bars (n=3), and thus more reproducible data, when compared to those grown manually (**Figure 6B**).




Figure 5. Drug response assay workflow.

Organoids from each source were resuspended at 90 organoids per  $\mu\text{l}$  in full growth medium supplemented with 5% hydrogel matrix. 25  $\mu\text{l}$  of this organoid suspension was dispensed into each well of a white, clear-bottomed 384 well plate (Greiner) pre-coated with 10  $\mu\text{l}$  matrix. Organoids from frozen stocks of bioreactor expanded organoids were given a 48 hour recovery period before treatment with compounds, while organoids harvested directly from manual culture were treated immediately. 11 point compound titrations (2 or 3 fold) were applied in quadruplicate alongside DMSO controls. After 5 days, cell viability assessments were made using an end-point Cell Titer Glo 3D assay. Experiments were performed in triplicate, from three independent bioreactor or manual expansions.




Figure 6. Drug response data.

BrC-001 organoid response to a variety of compounds was measured using Cell Titer Glo 3D on day 5 of treatment. Data was normalized to DMSO control conditions (100%) and analysed by non-linear regression analysis. Resulting dose-response curves were used to calculate EC50 values for each drug, for each organoid expansion condition (n=3). **(A)** Heatmap of calculated EC50 values for all drugs and organoid expansion conditions. Lower values are depicted by blue colour, while higher values are depicted by red colour. Further analysis of curve fitting was performed, with \* indicating pairs of responses that fit to the same curve, and are therefore not statistically different from one another. **(B)** Dose response curves following treatment with (i) AZD8055 or (ii) Paclitaxel, expressed as % response compared to control DMSO wells (mean  $\pm$ SD, n=3 for each condition)

## Conclusion

When taken together, the results of this study confirm that Cellesce's bioprocess is compatible with the large scale expansion of breast cancer organoid lines. Moreover, batches of organoids generated in this way produce more reproducible data than those produced by more labour intensive and costly manual methods.

Following the success of this study, Cellesce is currently working alongside Cardiff University to generate large scale batches of organoid lines from a variety of molecular subtypes of breast cancer, with the first lines available in late 2020.