



Hyperion XTi™ Imaging System

Uniquely designed for clinical research

Taking the complexity out of spatial biology

Introducing the **fastest and most reliable workflow** for high-plex imaging.

A next-generation solution using proven Imaging Mass Cytometry™ (IMC™) technology

- Faster detection for higher throughput
- Lower limit of detection to visualize dim markers
- Enhanced whole slide panorama to select region of interest in minutes
- Improved dual suspension and imaging modes available


HIGHLIGHTS

The only technology with

- 40-plus markers imaged simultaneously
- No autofluorescence interference
- Protein and RNA co-detection
- Integrated cell segmentation
- All-at-once batch slide staining and storage

From tissue collection to high-plex insights in 3 days*

Simplify your spatial imaging workflow without time-consuming panel design, acquisition cycles and management of spectral overlap.




Explore mechanism of action, disease progression and therapeutic outcomes **with confidence.**

*After panel and image analysis optimization

IMC applies **purified heavy metal labels**, not normally found in biological systems, instead of fluorophores.


40-plus markers without compromise

IMC is the only technology equipped to handle high-plex imaging for all tissue types – including lung, bone marrow and brain – without the interference of autofluorescence.



Capture targets of interest – in a single run


Simultaneously capture cell subpopulations, complex networks, cell-to-cell dynamics, cell distribution and biomarker expression.



Tissue architecture
Protein modifications
Signaling pathway activation
Cell injury states
Cell proliferation
Transcriptional signature
Cytokines and chemokines
And more ...


High-plex panel design – with ease

Start with ready-to-go panels and easily customize targets of interest. **Design your panel in days, not months.**



Established leader

For peer-reviewed high-plex imaging publications.



As of December 2022. A review of the Akoya Biosciences® website shows 12 peer-reviewed publications for PhenoCycler® (formerly CODEX®). The Ionpath® website lists 23 publications, but only 12 are peer-reviewed.

Sources: akoyabio.com/publication/; ionpath.com/publications/

Simplify spatial imaging with IMC
[standardbio.com/imc](https://www.standardbio.com/imc)

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Limited Use Label License:** The purchase of this Standard BioTools Instrument and/or Consumable conveys to the purchaser the limited, nontransferable right to use only with Standard BioTools Consumables and/or Instruments respectively except as approved in writing by Standard BioTools Inc. (f.k.a. Fluidigm Corporation): www.standardbio.com/legal/salesterms. **Patents:** www.standardbio.com/legal/notices. **Trademarks:** Standard BioTools, the Standard BioTools logo, Fluidigm, the Fluidigm logo, Hyperion, Hyperion XTi, Imaging Mass Cytometry, IMC and XTi are trademarks and/or registered trademarks of Standard BioTools Inc. or its affiliates in the United States and/or other countries. All other trademarks are the sole property of their respective owners. ©2023 Standard BioTools Inc. All rights reserved. 02/2023

FLDM-01089 Rev 01