

Explore the Power of More with high-plex protein profiling. More proteins, more pathways, more pQTLs, more proof.

Largest menu

7,000 or 11,000 protein measurements in a 55 µL biological sample*

Reproducibility

~5% median coefficient of variation (CV)

Sensitivity

10-log dynamic range from fM to uM

Specificity


Aptamer-based SOMAmer® Reagents designed for epitope binding

Publications

800+ and growing

www.somalogic.com/ publications/

Realize every sample's full potential: human plasma and serum Detect proteins across every major biological pathway



^{*}SomaScan 11K Assay: plasma and serum *SomaScan 7K Assay: plasma, serum, urine, CSF and other sample types

Identify biomarkers and elucidate pathways in discovery and translational science:

- Target identification and validation
- · Mechanism of action
- · Diesease biology modeling
- Pharmacokinetics/pharmacodyamics
- Enrichment of clinical trials

- · Stratification of clinical trials
- · Disease sub-phenotyping
- Physiology investigation
- Pathway regulation investigation
- · Dovetailing with clinical care

Key disease areas

Cardiovascular disease

Neurology

Diabetes

Nonalcoholic fatty liver disease

Oncology & Immunotherapy

Rare disease

SomaScan 7K Assay: Custom and disease specific panels and SomaSignal® tests

are also available for a more targeted approach to your protein profiling.

How to work with our CLIA-certified lab in the U.S.



-2.5 0.0 2.5

Initiate SomaScan Assay Services

- Gain access to our Life Science customer portal and initiate your project
- Choose from our SomaScan 11K Assay or the SomaScan 7K Assay which includes customized and disease focused panel
- Submit your samples to SomaLogic.
 Upon receipt, samples are processed and prepared for assaying at our CLIA/CAP-certified lab

The SomaScan Assay

- Before being processed, your samples are diluted and prepared with SomaScan Assay components including our engineered SOMAmer protein capture reagents
- Samples are transferred for hybridization to a custom microarray; data are read out upon completion of the SomaScan binding assay

Data Generation

- Data files are generated, formatted, and provided to you through our confidential customer portal
- Data can be visualized using tools in the SomaLogic portal or can be processed by your team

SomaLogic's project managers provide additional support through the process.

Abbreviations: CAP, College of American Pathologists; CLIA, Clinical Laboratory Improvement Amendments.

"SomaScan [Assay] is the state-of-art proteomic platform in multiplex, sensitivity and throughput."

Customer response to tech validate survey

Learn more about the SomaScan Assay somalogic.com/ somascan-assay-services/ ISO certifications
ISO 15189 Accreditation
Certificate
ISO 27001 – Certificate
Reissue – 2022
ISO 9001 Certificate

SL00000280 Rev 5 2024-01

SomaScan Assay Services Brochure



